

**ATCC medium: 2221 Quarter-strength cornmeal agar**

Cornmeal Infusion (see below)....250.0 ml  
Agar.....15.0 g  
Distilled water.....750.0 ml

Boil to dissolve agar. Autoclave at 121C for 15 minutes.

*Cornmeal Infusion:*

Yellow cornmeal.....50.0 g  
Distilled water.....1.0 L

Add cornmeal to the distilled water. Heat and bring to a boil. Simmer for 10 minutes. Filter through cheesecloth. Return volume to 1 liter.