

ATCC medium: 2235 PG medium for *Halosimplex*

NaCl.....200.0 g
MgSO₄ . 7H₂O20.0 g
KCl.....4.0 g
K₂HPO₄0.5 g
(NH₄)₂SO₄.....1.0 g
Sodium pyruvate*.....2.5 g
Glycerol.....2.5 g
Distilled water.....1.0 L

*Sodium acetate may be used instead of sodium pyruvate.
Adjust pH to 7.4. Autoclave at 121C for 15 minutes.