

ATCC medium: 2236 GM medium with 4-Chlorophenol

K ₂ HPO ₄	2.1 g
KH ₂ PO ₄	0.4 g
NH ₄ NO ₃	0.5 g
MgSO ₄ . 7H ₂ O	0.2 g
CaCl ₂ . 2H ₂ O	0.023 g
FeCl ₃ . 6H ₂ O solution (1 mg/ml)	2.0 ml
Yeast extract.....	5.0 g
4-Chlorophenol.....	0.15 g
Distilled water.....	1.0 L

Combine all ingredients except for FeCl₃ solution. Autoclave at 121C for 15 minutes. Add filter-sterilized FeCl₃ solution to cooled medium.