

ATCC Medium: 2241 YEPD w/ 200mcg/ml Geneticin

Complete Medium

*YPD Broth (BD 242820).....	50 g
Geneticin Solution (add aseptically, see below).....	50 ml
DI Water.....	950 ml

*This product is commercially available from BD. To make agar product, use 65 g of YPD Agar (BD 242720) in place of YPD Broth.

Geneticin Solution

Geneticin.....	0.2 g
DI Water.....	50 ml

Add to base medium aseptically.

Note: If commercial product is unavailable, make YPD from scratch, following the formulation below:

Yeast Extract.....	10.0 g
Peptone.....	20.0 g
Dextrose.....	20.0 g
Agar.....	15.0 g
DI Water.....	1000 ml

Combine ingredients, heat to dissolve agar if necessary, and autoclave at 121°C.