

ATCC Medium 0023: Nutrient Agar/Broth w/ 5% Urea

Agar Base Medium

Nutrient Agar.....	23.0 g
25% Urea Solution (See Below).....	200.0 mL
DI Water.....	800.0 mL

Broth Base Medium

Nutrient Broth.....	8.0 g
25% Urea Solution (See Below).....	200.0 mL
DI Water.....	800.0 mL

Adjust pH to **6.8 +/- 0.2**. Heat to dissolve agar/broth.

Autoclave at 121°C.

25% Urea Solution

Urea.....	50.0 g
DI Water.....	200.0 mL

Filter sterilize and add to sterile base medium.