

**ATCC Medium: 2854 Nutrient Agar/Broth with 300 mcg/ml Ampicillin**

**Agar Medium**

Nutrient Agar (BD 213000).....23 g  
DI Water.....950 ml

Autoclave at 121°C and allow to cool down to 45-50°C. Aseptically add the following filter sterilized solution:

**Ampicillin Solution**

Ampicillin.....300 mg  
DI Water.....50 ml

Dispense as required. **Store at 4°C. Do not store longer than 7-10 days.**

**Broth Medium**

Nutrient Broth (BD 234000).....8.0 g  
DI Water.....950 ml

Autoclave at 121°C and allow to cool down to 45-50°C. Aseptically add the following filter sterilized solution:

**Ampicillin Solution**

Ampicillin.....300 mg  
DI Water.....50 ml

Dispense as required. **Store at 4°C. Do not store longer than 7-10 days.**