

ATCC Medium: 2938 LB Agar/Broth, Miller w/200ug/ml Erythromycin

Agar Medium

LB Agar (BD 244520).....40 g
DI Water.....950 ml

Autoclave at 121°C.

Broth Medium

LB Broth (BD 244620).....25.0 g
DI Water.....950 ml

Autoclave at 121°C. Add following filter-sterilized solution once base media has cooled to 55°C.

Erythromycin Solution

Erythromycin.....200 mg
Ethanol.....5 ml
DI Water.....45 ml

Dissolve erythromycin in ethanol and then add water. Filter sterilize.