

ATCC medium: 304 Brooks agar

Cornmeal.....5.0 g
Distilled water to.....1.0 L

Simmer cornmeal in 500 ml water for 30 minutes; filter through cotton;
add water to filtrate to 1.0 L.

Add:

Yeast extract.....2.0 g
Casitone (BD 225930).....1.0 g
Malt extract.....1.0 g
 KH_2PO_4 1.5 g
 $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$ 0.5 g
Agar.....20.0 g

Autoclave at 121C for 15 minutes.