

**ATCC Medium: 336 Potato Dextrose Agar (PDA)**

Potato Dextrose Agar (BD 213400).....39.0 g  
DI Water.....1000 ml

Autoclave at 121°C.

\*This product is commercially available from BD.

To make medium from scratch, follow formulation below:

**Potato Dextrose Agar Composition**

Diced potatoes.....300.0 g  
Glucose.....20.0 g  
Agar.....15.0 g  
DI Water.....1000 ml

Boil finely diced potatoes in 500 ml of water until thoroughly cooked (about 1 hour); filter through cheesecloth and add water to filtrate to 1000 ml. Heat the filtrate to dissolve the agar. Add the glucose before sterilization. Autoclave at 121°C.