

ATCC Medium: 34 Neomycin Agar

Beef Extract.....	1.5 g
Yeast Extract.....	3.0 g
Peptone.....	6.0 g
Glucose.....	1.0 g
Casitone.....	4.0 g
Neomycin Solution (See Below).....	50 ml
Agar.....	15.0 g
DI water.....	950 ml

Adjust to pH 7.0.
Autoclave at 121°C.

Neomycin Solution

Neomycin.....	1.0 g
DI Water.....	50 ml

*Filter sterilize solution and add to sterile base medium.