

ATCC medium: 405 BSL for *Corynebacterium*

NH₄Cl1.0 g
MgSO₄ . 7H₂O0.2 g
KH₂PO₄3.0 g
Na₂HPO₄6.0 g
Soytone (BD 243620).....10.0 g
Agar (if needed).....15.0 g
Distilled water.....1.0 L

Adjust pH to 6.8-7.2. Separately autoclave a 10% (w/v) solution of lactose and aseptically add to above for a final concentration of 1% lactose.

Autoclave at 121C for 15 minutes.