

ATCC Medium: 415 Potato Extract Agar/Broth

For Solid Media:

Nutrient Agar23.0 g
Potato Extract (see below).....20.0 ml
DI Water.....980.0 mL

Autoclave at 121°C

To prepare potato extract, boil 300 g peeled, diced potatoes in 500 ml distilled water for half an hour. Decant the infusion and use in the above medium. (Put through cheese cloth a few times.)

For Broth

Nutrient Broth.....8.0 g
Potato Extract (see above)20.0 mL
DI Water.....980.0 mL

Autoclave at 121°C