

**ATCC medium: 490 Brucella broth with additives**

ATCC Medium 489 (see below) with 30 g NaCl per liter of medium.

*ATCC Medium 489:*

Brucella Broth (BD 211088).....23.8 g  
Agar, Noble (BD 214230) (if needed).....15.0 g  
Distilled water.....850.0 ml

Autoclave at 121C for 15 minutes. Cool to 50C and aseptically add the following filter-sterilized components:

Horse serum (heat-inactivated).....100.0 ml  
Yeast Extract Solution (GIBCO 18180).....50.0 ml