

ATCC medium: 608 Na/suc medium with Penicillin G

Basal medium:

Brain Heart Infusion (BD 237500)	37.0 g
Yeast extract	5.0 g
NaCl	9.3 g
Sucrose	97.3 g
MgSO ₄	0.25 g
Agar (if needed)	13.3 g
Deionized water to	1.0 L

Adjust pH to 7.4. Autoclave at 121C for 15 minutes. To 900 ml of this solution aseptically add:

Horse serum, gamma globulin-free (heat-inactivated)	100.0 ml
Penicillin G (Pfizer)	50,000 U