

**ATCC Medium 0661: Alkaline Bacillus Medium**

**Base Media:**

|                                       |          |
|---------------------------------------|----------|
| Peptone.....                          | 10.0 g   |
| Yeast Extract.....                    | 5.0 g    |
| Glucose.....                          | 10.0 g   |
| K <sub>2</sub> HPO <sub>4</sub> ..... | 1.0 g    |
| Agar.....                             | 15.0 g   |
| DI Water.....                         | 900.0 mL |

*Sterilize at 121°C for 15 minutes.*

*Add the following solution to the base media:*

**10% Na<sub>2</sub>CO<sub>3</sub> Solution:**

|                                       |          |
|---------------------------------------|----------|
| Na <sub>2</sub> CO <sub>3</sub> ..... | 10.0 g   |
| DI Water.....                         | 100.0 mL |

*Filter sterilize.*