

ATCC medium: 717 PFS medium

Peptone, Bacto (BD 211677).....10.0 g
Fumaric acid.....2.0 g
(NH₄)₂SO₄.....1.0 g
MnSO₄ . H₂O0.2 mg
MgSO₄ . 7H₂O0.5 g
FeCl₃ . 6H₂O0.2 mg
Distilled water.....1.0 L

Adjust pH to 7.0 with KOH. Incorporate 0.15% agar for maintenance medium. Autoclave at 15 lbs for 15-20 minutes.