

ATCC medium: 731 Luminous medium

NaCl.....	30.0 g
K ₂ HPO ₄	3.9 g
KH ₂ PO ₄	2.1 g
NH ₄ Cl	5.0 g
Yeast extract.....	5.0 g
Tryptone.....	5.0 g
1 M Tris buffer, pH 7.5.....	50.0 ml
Glycerol.....	3.0 ml
MgSO ₄ . 7H ₂ O	1.0 g
KCl.....	0.75 g
CaCO ₃	1.0 g
Agar (if needed).....	20.0 g
Distilled water to.....	1.0 L

Adjust pH to 7.2.

Autoclave 15 minutes at 121C.