

ATCC medium: 733 SCY medium

Solution A:

Trypticase Soy Broth without Dextrose (BD 211774).....0.25 g
Casitone (BD 225930).....0.75 g
Agar.....10.0 g
Distilled water.....1.0 L

Adjust pH to 7.0. Autoclave 15 minutes at 121C. Slant tubes.

Solution B:

Sucrose.....1.0 g
Vitamin B12.....0.01 mg
Thiamine.....0.4 mg
Yeast Extract (BD 212750).....0.25 g
Tap water (pH 8.5).....100.0 ml

Filter-sterilize. Pipette 1 to 2 ml of Solution B onto slants of Solution A after inoculating culture.