

ATCC medium: 743 Alkvisco medium

Beef Extract (BD 212610).....10.0 g
Peptone.....10.0 g
NaCl.....5.0 g
KCN.....10.0 mg
Acrylonitrile*.....0.5 g
Agar.....15.0 g
Distilled water.....1.0 L

Adjust pH to 6.5-8.0.

*Filter-sterilize the acrylonitrile solution and add to the sterile basal medium which has been autoclaved at 121C for 10 minutes.

Acrylonitrile is a carcinogen; use appropriate precautions.