

ATCC medium: 782 Synthetic seawater medium

NaCl.....27.0 g
CaCl₂.....0.3 g
KCl.....0.6 g
MgSO₄ . 7H₂O7.0 g
Sodium glycerophosphate.....0.2 g
Tris buffer.....2.0 g
Monosodium glutamate.....5.0 g
Glucose.....1.0 g
Vitamin B12.....1.0 mcg
Distilled water.....1.0 L

Adjust pH to 7.5 with N HCl.
Autoclave at 121C for 15 minutes.