

ATCC medium: 805 CM plus ye medium

Casamino Acids, Vitamin Assay (BD 228830).....	7.5 g
Yeast Extract (BD 212750).....	10.0 g
MgSO ₄ . 7H ₂ O	20.0 g
Sodium citrate.....	3.0 g
KCl.....	2.0 g
NaCl.....	200.0 g
Fe*.....	10.0 ppm
Distilled water.....	1.0 L

Adjust pH to 7.4 with NaOH and autoclave at 15 lbs for 10 minutes. For a solid medium add 1.5% agar.

*Prepare 4.98% FeSO₄ . 7H₂O solution; add 1 ml per liter of medium.