

ATCC medium: 822 C-3-N *Spiroplasma* medium

PPLO Broth w/o CV (BD 255420).....1.5 g
Sucrose.....12.0 g
Phenol red.....0.001 %
Distilled water.....69.0 ml

Autoclave at 15 lbs. for 15 minutes. Add aseptically to above:

Horse serum.....20.0 ml
Yeast Extract Solution (GIBCO 18180).....5.0 ml
CMRL Medium (Flow Labs. 2-014D).....0.5 ml

Adjust pH to 7.5.