

**ATCC medium: 907 Cellulose agar**

Ammonium sulfate.....0.5 g  
L-Asparagine.....0.5 g  
KH<sub>2</sub>PO<sub>4</sub> .....1.0 g  
KCl.....0.5 g  
MgSO<sub>4</sub>.....0.2 g  
CaCl<sub>2</sub>.....0.1 g  
Yeast extract.....0.5 g  
Agar.....20.0 g  
Cellulose (ball milled).....10.0 g  
Distilled water up to.....1.0 L

Autoclave at 121C for 15 minutes.  
Adjust pH to 6.2.