

ATCC medium: 930 Medium N

Distilled water.....	1.0 L
Yeast Nitrogen Base w/o Amino Acids (BD 291940).....	6.7 g
Casamino Acids, Vitamin Assay (BD 228830).....	2.0 g
Glucose.....	20.0 g
Agar.....	20.0 g
Isoleucine.....	0.1 g
Valine.....	0.1 g

Autoclave at 121C for 15 minutes.

After the autoclaved Medium N has been allowed to cool down to 60C, aseptically add 5'-dTMP filter-sterilized at the rate of 15 mcg/ml.