

ATCC Medium: 96 Nitrobacter Medium B

NaNO ₂	1.0 g
K ₂ HPO ₄	0.5 g
NaCl.....	0.3 g
MgSO ₄	0.5 g
MnSO ₄	0.002 g
Fe ₂ (SO ₄) ₃	0.005 g
DI Water.....	1000 ml

Adjust pH to 7.5. Dispense in 100 ml amounts in flasks. Sterilize. Wash marble chips (Fisher M123-3) in distilled water. Place several chips in 18 x 150 mm test tubes. Sterilize one hour. Aseptically add medium to cover the chips.