

ATCC Medium 9: Nutrient Agar/Broth w/ 2% Yeast Extract

Nutrient Agar (BD 213000)	23.0 g
Nutrient Broth (BD 234000).....	8.0 g
Yeast Extract.....	20.0 g
DI Water.....	1000.0 mL

Autoclave at 121°C.

Dispense into appropriate vessel type.

To make medium from scratch, follow formulation below:

BD Nutrient Agar/Broth Composition

Beef Extract.....	3.0 g
Peptone.....	5.0 g
*Agar.....	15.0 g

*Omit agar for broth medium

Final pH 6.8 +/- 0.2.