

ATCC Medium: 1525 Seawater Sonneborn's Paramecium Medium

2X Seawater (Seven Seas).....	500.0	ml
*2X Modified ATCC Medium 802.....	500.0	ml

Filter-sterilize seawater. Add seawater to cooled autoclaved 2X Modified ATCC Medium 802.

***2X Modified ATCC Medium 802**

Rye Cerophyl.....	5.0	g
Distilled water.....	1000.0	ml

Add rye cerophyl to distilled water and boil for 5 minutes. Add distilled water to compensate for evaporation. Filter through Whatman #1 filter paper. Autoclave at 121 C.