

ATCC medium: 1620 CH-1 medium

Yeast extract.....5.0 g
Hy-Case SF.....5.0 g
NaCl.....250.0 g
Tris.....12.0 g
Solution 1 (see below).....50.0 ml
Glycerol.....10.0 g
Distilled water to.....1.0 L

Adjust for a final pH of 7.4. Autoclave at 121C for 15 minutes.

Solution 1:

MgCl₂ . 6H₂O40.0 g
KCl.....4.0 g
CaCl₂ . 2H₂O0.4 g
Distilled water.....1.0 L