

ATCC medium: 1634 LB medium (ATCC medium 1065) with 50 mcg/ml kanamycin

ATCC Medium 1065 (see below) with 50.0 mcg/ml kanamycin

Prepare a filter-sterilized solution of kanamycin and aseptically add enough of it to the sterile basal medium to achieve the stated concentration.

ATCC Medium 1065:

Tryptone (BD 211705).....10.0 g
Yeast Extract (BD 212750)....5.0 g
NaCl.....10.0 g
Distilled water.....1.0 L

Autoclave at 121C for 15 minutes.