

ATCC Medium 1685: NYC III

Base Medium:

HEPES.....	4.0 g
Proteose Peptone No. 3.....	15.0 g
NaCl	5.0 g
Distilled water	875.0 mL

Adjust for final pH 7.3.

Autoclave at 121° C for 15 minutes.

Cool to room temperature and aseptically add the following additive solution.

Additive Solution:

Glucose	5.0 g
Fresh yeast extract (Gibco 360-8180)	25.0 mL
Heat Inactivated Horse Serum	100.0 mL

Mix well and filter sterilize.