**Product Sheet** 

# Genomic DNA from Salmonella enterica subsp. enterica serovar Thompson strain MZ1490

**BAA-1738D-5<sup>™</sup>** 

### Description

Genomic DNA isolated from Salmonella enterica subsp. enterica serovar Thompson strain MZ1490. This bacterial strain is also available as ATCC<sup>®</sup> Catalog No. BAA-1738<sup>™</sup>. **Organism:** Salmonella enterica **Derived from:** Salmonella enterica subsp. enterica serovar Thompson MZ1490 [SARB62] (ATCC BAA-1738) **Genome sequenced strain:** Yes **Type strain:** No **Mass:** 5 µg **Shipping information:** Stored in 1X TE buffer

Storage Conditions Product format: Freeze-dried

### Intended Use

This product is intended for laboratory research use only. It is not intended for any animal or human therapeutic use, any human or animal consumption, or any diagnostic use.

BSL 1



# Genomic DNA from Salmonella enterica subsp. enterica serovar Thompson strain MZ1490

### BAA-1738D-5

ATCC determines the biosafety level of a material based on our risk assessment as guided by the current edition of *Biosafety in Microbiological and Biomedical Laboratories (BMBL)*, U.S. Department of Health and Human Services. It is your responsibility to understand the hazards associated with the material per your organization's policies and procedures as well as any other applicable regulations as enforced by your local or national agencies.

### **Certificate of Analysis**

For batch-specific test results, refer to the applicable certificate of analysis that can be found at www.atcc.org.

### Handling Procedures

Centrifuge tube prior to opening to prevent loss of pelleted material

- 1. Rehydrate contents of vial with molecular grade  $H_2O$ .
- 2. Place vial at 37°C for 1 hour or at 2°C to 8°C overnight.
- 3. For more complete rehydration and to fully recover DNA, incubate the sample overnight at 4°C while rocking; then incubate for 1 hour at 65°C. Resuspending the dried DNA in  $\ge$  250 µL may give better results.

### **Quality Control Specifications**

**Total amount:** Total DNA by PicoGreen<sup>®</sup> measurement was found to be approximately 5  $\mu$ g.

**Integrity:** Integrity of DNA was determined by electrophoresis on a 1% agarose gel stained with SYBR Safe<sup>™</sup>, and was found to be of high molecular weight.

**Functional tests:** Functional activity was confirmed by PCR amplification of the 16S ribosomal RNA gene.

**Identity:** Identity confirmed by sequencing of 16S ribosomal RNA gene (first ~500 base pairs).



### Genomic DNA from Salmonella enterica subsp. enterica serovar Thompson strain MZ1490 BAA-1738D-5

Notes

Genomic DNA is appropriate for PCR and other molecular biology applications.

### **Material Citation**

If use of this material results in a scientific publication, please cite the material in the following manner: Genomic DNA from *Salmonella enterica* subsp. *enterica* serovar Thompson strain MZ1490 (ATCC BAA-1738D-5)

### References

References and other information relating to this material are available at www.atcc.org.

### Warranty

The product is provided 'AS IS' and the viability of ATCC<sup>®</sup> products is warranted for 30 days from the date of shipment, provided that the customer has stored and handled the product according to the information included on the product information sheet, website, and Certificate of Analysis. For living cultures, ATCC lists the media formulation and reagents that have been found to be effective for the product. While other unspecified media and reagents may also produce satisfactory results, a change in the ATCC and/or depositor-recommended protocols may affect the recovery, growth, and/or function of the product. If an alternative medium formulation or reagent is used, the ATCC warranty for viability is no longer valid. Except as expressly set forth herein, no other warranties of any kind are provided, express or implied, including, but not limited to, any implied warranties of merchantability, fitness for a particular purpose, manufacture according to cGMP standards, typicality, safety, accuracy, and/or noninfringement.



### Genomic DNA from Salmonella enterica subsp. enterica serovar Thompson strain MZ1490 BAA-1738D-5 Disclaimers

This product is intended for laboratory research use only. It is not intended for any animal or human therapeutic use, any human or animal consumption, or any diagnostic use. Any proposed commercial use is prohibited without a license from ATCC.

While ATCC uses reasonable efforts to include accurate and up-to-date information on this product sheet, ATCC makes no warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. ATCC does not warrant that such information has been confirmed to be accurate or complete and the customer bears the sole responsibility of confirming the accuracy and completeness of any such information.

This product is sent on the condition that the customer is responsible for and assumes all risk and responsibility in connection with the receipt, handling, storage, disposal, and use of the ATCC product including without limitation taking all appropriate safety and handling precautions to minimize health or environmental risk. As a condition of receiving the material, the customer agrees that any activity undertaken with the ATCC product and any progeny or modifications will be conducted in compliance with all applicable laws, regulations, and guidelines. This product is provided 'AS IS' with no representations or warranties whatsoever except as expressly set forth herein and in no event shall ATCC, its parents, subsidiaries, directors, officers, agents, employees, assigns, successors, and affiliates be liable for indirect, special, incidental, or consequential damages of any kind in connection with or arising out of the customer's use of the product. While reasonable effort is made to ensure authenticity and reliability of materials on deposit, ATCC is not liable for damages arising from the misidentification or misrepresentation of such materials.

Please see the material transfer agreement (MTA) for further details regarding the use of this product. The MTA is available at www.atcc.org.

### Copyright and Trademark Information

© ATCC 2023. All rights reserved.

ATCC is a registered trademark of the American Type Culture Collection.



# Genomic DNA from Salmonella enterica subsp. enterica serovar Thompson strain MZ1490

BAA-1738D-5 **Revision** 

This information on this document was last updated on 2024-10-25

## **Contact Information**

ATCC 10801 University Boulevard Manassas, VA 20110-2209 USA US telephone: 800-638-6597 Worldwide telephone: +1-703-365-2700 Email: tech@atcc.org or contact your local distributor

