



Genomic DNA from *Porphyromonas gingivalis* strain W83

BAA-308D-5™

Description

Genomic DNA isolated from *Porphyromonas gingivalis* strain W83. This bacterial strain is also available as ATCC® BAA-308™.

Organism: *Porphyromonas gingivalis* (Coykendall et al.) Shah and Collins

Derived from: *Porphyromonas gingivalis* W83 (ATCC BAA-308)

Genome sequenced strain: Yes

Type strain: No

Mass: 5 µg

Shipping information: Stored in 1X TE buffer

Storage Conditions

Product format: Dried

Intended Use

This product is intended for laboratory research use only. It is not intended for any animal or human therapeutic use, any human or animal consumption, or any diagnostic use.

BSL 1

ATCC determines the biosafety level of a material based on our risk assessment as guided by the current edition of *Biosafety in Microbiological and Biomedical Laboratories (BMBL)*, U.S. Department of Health and Human Services. It is your responsibility to

Genomic DNA from *Porphyromonas gingivalis* strain W83

Product Sheet

BAA-308D-5

understand the hazards associated with the material per your organization's policies and procedures as well as any other applicable regulations as enforced by your local or national agencies.

Certificate of Analysis

For batch-specific test results, refer to the applicable certificate of analysis that can be found at www.atcc.org.

Handling Procedures

Centrifuge tube prior to opening to prevent loss of pelleted material

1. Rehydrate contents of vial with molecular grade H₂O.
 2. Place vial at 37°C for 1 hour or at 2°C to 8°C overnight.
 3. For more complete rehydration and to fully recover DNA, incubate the sample overnight at 4°C while rocking; then incubate for 1 hour at 65°C. Resuspending the dried DNA in ≥ 250 µL may give better results.
-

Quality Control Specifications

Integrity: Integrity of DNA was determined by electrophoresis on a 1% agarose gel stained with ethidium bromide, and was found to be intact and of high molecular weight.

Notes

Genomic DNA isolated from bacteria is appropriate for PCR and other molecular biology applications.

Material Citation

Genomic DNA from *Porphyromonas gingivalis* strain W83

BAA-308D-5

If use of this material results in a scientific publication, please cite the material in the following manner: Genomic DNA from *Porphyromonas gingivalis* strain W83 (ATCC BAA-308D-5)

References

References and other information relating to this material are available at www.atcc.org.

Warranty

The product is provided 'AS IS' and the viability of ATCC® products is warranted for 30 days from the date of shipment, provided that the customer has stored and handled the product according to the information included on the product information sheet, website, and Certificate of Analysis. For living cultures, ATCC lists the media formulation and reagents that have been found to be effective for the product. While other unspecified media and reagents may also produce satisfactory results, a change in the ATCC and/or depositor-recommended protocols may affect the recovery, growth, and/or function of the product. If an alternative medium formulation or reagent is used, the ATCC warranty for viability is no longer valid. Except as expressly set forth herein, no other warranties of any kind are provided, express or implied, including, but not limited to, any implied warranties of merchantability, fitness for a particular purpose, manufacture according to cGMP standards, typicality, safety, accuracy, and/or noninfringement.

Disclaimers

This product is intended for laboratory research use only. It is not intended for any animal or human therapeutic use, any human or animal consumption, or any diagnostic use. Any proposed commercial use is prohibited without a [license from ATCC](#).

Genomic DNA from *Porphyromonas gingivalis* strain W83

Product Sheet

BAA-308D-5

While ATCC uses reasonable efforts to include accurate and up-to-date information on this product sheet, ATCC makes no warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. ATCC does not warrant that such information has been confirmed to be accurate or complete and the customer bears the sole responsibility of confirming the accuracy and completeness of any such information.

This product is sent on the condition that the customer is responsible for and assumes all risk and responsibility in connection with the receipt, handling, storage, disposal, and use of the ATCC product including without limitation taking all appropriate safety and handling precautions to minimize health or environmental risk. As a condition of receiving the material, the customer agrees that any activity undertaken with the ATCC product and any progeny or modifications will be conducted in compliance with all applicable laws, regulations, and guidelines. This product is provided 'AS IS' with no representations or warranties whatsoever except as expressly set forth herein and in no event shall ATCC, its parents, subsidiaries, directors, officers, agents, employees, assigns, successors, and affiliates be liable for indirect, special, incidental, or consequential damages of any kind in connection with or arising out of the customer's use of the product. While reasonable effort is made to ensure authenticity and reliability of materials on deposit, ATCC is not liable for damages arising from the misidentification or misrepresentation of such materials.

Please see the material transfer agreement (MTA) for further details regarding the use of this product. The MTA is available at www.atcc.org.

Copyright and Trademark Information

© ATCC 2023. All rights reserved.

ATCC is a registered trademark of the American Type Culture Collection.

Revision

This information on this document was last updated on 2024-10-25

Genomic DNA from *Porphyromonas gingivalis* strain W83

BAA-308D-5

Contact Information

ATCC

10801 University Boulevard

Manassas, VA 20110-2209

USA

US telephone: 800-638-6597

Worldwide telephone: +1-703-365-2700

Email: tech@atcc.org or contact your local distributor
