

**ATCC Medium: 1395 Middlebrook 7H9 Broth w/ADC Enrichment**

Middlebrook 7H9 Broth.....	4.7g
Glycerol.....	2.0 ml
DI Water.....	900 ml

Mix well. Autoclave at 121°C for 10 minutes.

Aseptically add 100ml ADC Enrichment:

**ADC Enrichment**

Bovine Serum Albumin, Frac V.....	5.0 g
Dextrose.....	2.0 g
Catalase (beef).....	3 mg
DI Water.....	100ml

Filter Sterilize.

Store any unused portion with lot information and Media Number in 2-8°C.