

**ATCC medium: 1176 *Halobacterium* medium**

NaCl.....156.0 g  
MgCl<sub>2</sub> . 6H<sub>2</sub>O .....13.0 g  
MgSO<sub>4</sub> . 7H<sub>2</sub>O .....20.0 g  
CaCl<sub>2</sub> . 6H<sub>2</sub>O .....1.0 g  
KCl.....4.0 g  
NaHCO<sub>3</sub>.....0.2 g  
NaBr.....0.5 g  
Yeast extract.....5.0 g  
Glucose.....1.0 g  
Distilled water to.....1.0 L

Adjust pH to 7.0. Autoclave at 121C for 15 minutes.