

ATCC Medium 1063: Yeast-Tryptone Medium with 200mg/L of Streptomycin

Base Medium

Yeast Extract	5.0 g
Tryptone	10.0 g
NaCl	10.0 g
Streptomycin solution (see below)	200.0 mg
DI Water	980.0 mL

Autoclave base medium at 121°C for 15 minutes.

Aseptically add streptomycin solution.

Dispense into appropriate vessel type.

Streptomycin solution

Streptomycin	200.0 mg
DI Water	20.0 mL

Filter sterilize.