

ATCC medium: 2622 Halophilic medium

Yeast extract.....	0.10 g
MgCl ₂ . 6H ₂ O.....	20.00 g
K ₂ SO ₄	5.00 g
CaCl ₂ . 2H ₂ O.....	0.10 g
NH ₄ Cl.....	0.50 g
KH ₂ PO ₄	0.05 g
NaCl.....	180.00 g
Sucrose.....	0.50 g
HEPES.....	6.00 g
Agar.....	20.00 g
Distilled water.....	1.0 L

Adjust the pH to 6.5. Autoclave at 121°C for 15 minutes.