

**101**      **Photobacterium Broth**

Difco Photobacterium Broth has been discontinued.

Tryptone	5.0	g
Yeast Extract	2.5	g
NH <sub>4</sub> Cl	0.3	g
MgSO <sub>4</sub> x7H <sub>2</sub> O	0.3	g
FeCl <sub>3</sub> x6H <sub>2</sub> O	0.01	g
CaCO <sub>3</sub>	1.0	g
KH <sub>2</sub> PO <sub>4</sub>	3.0	g
Sodium Glycerophosphate	23.5	g
NaCl	30.0	g
Agar (if required)	15.0	g
DI water	1000	ml

Heat to dissolve completely. Boil if agar is included.

Adjust pH to 7.0 +/- 0.2.

Autoclave at 121 C for 15 minutes.