

ATCC medium: 329 Mineral salts agar

NaNO₃2.0 g
MgSO₄0.5 g
KCl.....0.5 g
Fe₂(SO₄)₃ . H₂O0.01 g
KH₂PO₄0.14 g
K₂HPO₄1.2 g
Agar.....15.0 g
Yeast extract.....0.02 g
Distilled water.....1.0 L

Autoclave at 121C for 15 minutes.
Adjust pH to 7.2. Place strip of sterile filter paper on cooled agar
slant and inoculate on filter strip.