

**ATCC Medium: 1638 LB medium (ATCC medium 1065) with 5 mcg/ml kanamycin**

ATCC Medium 1065 (see below) with 5.0 mcg/ml kanamycin

Prepare a filter-sterilized solution of kanamycin and aseptically add enough of it to the sterile basal medium to achieve the stated concentration.

*ATCC Medium 1065:*

Tryptone (BD 211705).....	10.0 g
Yeast Extract (BD 212750).....	5.0 g
NaCl.....	10.0 g
Distilled water.....	1.0 L

Autoclave at 121C for 15 minutes.