

ATCC medium: 807 Brain heart infusion blood agar

Agar:

Brain Heart Infusion (BD 237500).....37.0 g
Agar.....18.0 g
Distilled water.....1.0 L

Heat to boiling to dissolve the medium completely. Dispense 5.0 ml per 16 X 125 screw-capped test tube, autoclave for 25 minutes at 121C. Cool to 48C. Add 0.5 ml per tube of sterile defibrinated rabbit blood and slant. After slants have set, cover with 0.5 ml of Brain Heart Infusion (Difco 0037) broth.

Overlay:

Brain Heart Infusion (BD 237500).....37.0 g
Distilled water.....1.0 L