

ATCC Medium: 2167 Haemophilus Test Medium

Solution A

Mueller-Hinton Broth.....22.0 g
Yeast Extract.....5.0 g
DI Water.....900 ml

Autoclave at 121C.

Solution B

Hemin (dissolve in NaOH).....15.0 mg
Nicotinamide adenine dinucleotide.....15.0 mg
DI Water.....100 ml

Filter sterilize.

Complete medium: Add solutions A and B together. Adjust for final pH 7.3 +/- 0.1. Dispense as required.