

**ATCC medium: 484 PPLO Broth w/ Additives for Mycoplasma**

PPLO broth w/o CV (BD 255420).....	14.7 g
Arginine HCl.....	2.1 g
Agar, Noble (BD 214230) (if needed).....	10.0 g
DI Water.....	700.0 ml

Autoclave at 121C for 15 minutes. Cool to 50C and aseptically add the following filter-sterilized components:

Glutamine.....	1.46 g
Phenol red.....	20.0 mg
100X MEM Vitamins (GIBCO 11120-052).....	1.0 ml
Horse serum.....	200.0 ml
Yeast Extract Solution (GIBCO 18180).....	100.0 ml

Adjust medium, if necessary, for final pH 7.1.