

**ATCC Medium: 1970 V8 Rye Agar**

CaCO<sub>3</sub> .....0.2 g  
V8 Juice.....50.0 ml  
Agar.....20.0 g  
Rye broth (see below).....1000 ml

Autoclave at 121°C.

**Rye broth**

Whole rye grains.....50.0 g  
DI Water.....1100 ml

Soak grains in distilled water for 24-36 hours at 24°C. Autoclave at 121°C for 30 minutes. Filter through 4 layers of cheesecloth and adjust the final filtrate volume to 1000 ml.