

ATCC medium: 313 Czapek's agar

ATCC Medium 312 (see below) plus 170 g sucrose

ATCC Medium 312:

NaNO ₃	3.0 g
K ₂ HPO ₄	1.0 g
MgSO ₄ . 7H ₂ O	0.5 g
KCl	0.5 g
FeSO ₄ . 7H ₂ O	0.01 g
Sucrose (commercial grade)*	30.0 g
Agar	15.0 g
Distilled water	1.0 L

Combine all ingredients but sucrose; autoclave at 121C for 25 minutes.

*Add sucrose prior to dispensing.