

**ATCC medium: 485 PPLO broth or agar at pH 7.6 with additives for
*Mycoplasma***

PPLO Broth w/o CV (BD 255420).....14.7 g
Agar, Noble (BD 214230).....10.0 g
Distilled water.....700.0 ml

Autoclave at 121C for 15 minutes. Cool to 50C and aseptically add the
following filter-sterilized components:

L-Cystine . HCl . H₂O.....1.0 g
Horse Serum (not inactivated).....200.0 ml
Fresh baker's yeast extract (GIBCO 18180).....100.0 ml

Adjust medium, if necessary, for final pH 7.6.