

ATCC Medium: 0032 Green Top Agar

Yeast Extract.....	2.0 g
Tryptone.....	1.0 g
Sodium acetate.....	1.0 g
Soil Extract (See Below).....	50.0 ml
Agar.....	2.0 g
DI Water.....	950 ml

Adjust to pH 7.4 +/- 0.2

For solid medium, increase agar concentration to 15 grams/liter.

Autoclave at 121°C

Soil Extract

African Violet Soil.....	77.0 g
Na ₂ CO ₃	0.2 g
Tap Water.....	200 ml

Boil for 1 hour. Filter to clarify before use.