

ATCC medium: 1659 YA halophile medium

Na₂HPO₄3.8 g
KH₂PO₄1.3 g
(NH₄)₂SO₄.....1.0 g
NaCl.....100.0 g
Sodium acetate . 3H₂O10.0 g
Yeast extract.....1.0 g
Mg(NO₃)₂ . 6H₂O*.....1.0 g
Agar.....15.0 g
Distilled water to.....1.0 L

Adjust medium for final pH 7.2 with KOH. Autoclave at 121C for 15 minutes.

*Add aseptically after sterilization of medium.