

**ATCC Medium: 2758 Oatmeal Agar w/ KOH**

Oatmeal ("Quaker old fashioned" brand, large flakes).....	20.0 g
2N KOH.....	5.0 ml
Agar.....	20.0 g
DI Water.....	955 ml

Boil 20 g whole flake oatmeal in 500 ml water for 15 minutes. Strain through 4 layers of cheesecloth. Add 5 ml of 2N KOH. Bring volume to 1 L with distilled water. Add agar and boil. Autoclave at 121°C.