

ATCC medium: 757 YB medium

Methanol (filter-sterilized).....	20.0 ml
NaCl.....	3.0 g
Thiourea.....	0.1 g
Peptone.....	10.0 g
Yeast extract.....	5.0 g
Beef extract.....	7.0 g
Agar.....	20.0 g
Deionized water.....	980.0 ml

Combine all ingredients except methanol. Adjust pH to 7.2. Autoclave at 121C for 15 minutes.
Aseptically add filter-sterilized methanol to sterile basal medium.