

ATCC medium: 776 Mineral salts medium with methanol and yeast extract

NaNH ₄ HPO ₄ . 4H ₂ O.....	1.74 g
NaH ₂ PO ₄ . H ₂ O.....	0.54 g
KCl.....	0.04 g
MgSO ₄ . 7H ₂ O	0.2 g
FeSO ₄ . 7H ₂ O.....	5.0 mg
Yeast extract	0.2 g
Methanol.....	5.0 ml
Trace Minerals (see below).....	1.0 ml
Distilled water to.....	1.0 L

Combine all ingredients other than methanol. Adjust pH to 7.2.
Autoclave at 121C for 15 minutes. Add filter-sterilized methanol aseptically.

Trace Minerals:

H ₃ BO ₃	2.86 g
MnCl ₂ . 4H ₂ O	1.81 g
Na ₂ MoO ₄ . 2H ₂ O.....	25.0 mg
ZnSO ₄ . 7H ₂ O	0.22 g
CuSO ₄ . 5H ₂ O	0.08 g
CoCl ₂ . 6H ₂ O.....	0.06 g
Distilled water.....	1.0 L